

# Surginet

Elastic Bandage Net



## Appliaction

Burn , Sensitive Dermatopathy, Ventilation necessity in affected part



## Feature & Benefit

### ○ Economical and Efficient

- Cover anywhere with easy
- Unwoven strand
- Free cutting for purpose
- Washable and use repeatedly
- Consistent elasticity and no transformation

### ○ Convenient to Use

- Very simple to use especially for Head, Breast Shoulder, Waist and thigh
- Fix the gauge placed on the affected part and Patient can move easily

### ○ Enhance Remedial Value of Patient care.

- Ventilates the affected part
- Patient feels comfortable
- Prevents from side effects like dermatitis

## How to Use



1. Finger



2. Plam



3. Arm



4. Elbow



5. Knee



6. Foot



7. Thigh



8. head



9. Breast



6. Shoulder

# Surginet

## Surginet Standard Size

Section	Size	Use
0,5	1,1cm x 25m	More than two finger, Child's arm
1	1,6cm x 25m	More than two finger, Child's leg
2	2,0cm x 25m	Hand, Arm, Leg
3	2,2cm x 25m	Foot, Knee
4	2,5cm x 25m	Knee, Thigh, Child's head
5	3,0cm x 25m	Shoulder, thigh, Child's head
6	5,5cm x 25m	Breast, Waist, Body, Head
7	6,0cm x 25m	Breast, Waist, Body
8	6,5cm x 25m	Breast, Waist, Body



## Surgicap

### Appliaction

- Dressing for Head

### Feature & Benefit

- Easy to use
- Comfortable
- High vented material

### Size

- Large ( for adult )
- Small ( for children )



**WON INDUSTRY CO.**

222-5, Daeya-dong, Siheung-si, Gyeonggi, KOREA

TEL: (031)311-4024 FAX: (031)318-5573

E-mail: onesp77@naver.com

Surginet